

CryoStore™

Conical Freezing Bag

The OriGen CryoStore™ bag is indicated for freezing blood components.

Traceability

A permanent lot number is heat welded into the tail of each bag.

Durability

Better Material, Proven Performance

The hermetically sealed spike ports and film used in every CryoStore™ bag are made from a rugged EVA. Our customers have had great success with the strength of the bag over time.

Ease of Use

Designed In

OriGen can provide your choice of several tube sets (see reverse). Start a development project with us to reduce processing time and the possibility of contamination in the lab.



Technical Information

| Bag Size | Width, (cm) | Length, (cm) | Freeze volume min-max, ml | Fill volume (max), ml |
|----------|-------------|--------------|---------------------------|-----------------------|
| CSD400 | 12.7 | 20.3 | 60-135 | 420 |
| CSD2500 | 27.9 | 33.7 | 730-810 | 3000 |

CryoStore™

Conical Freezing Bag

General Information


Individually packed in Tyvek® pouch, supplied sterile, single use only.

Ordering Information:

For availability in your country, please contact your local distributor.


| Product Code | Description | Stock or Non-Stock | Case Qty |
|--------------|--|--------------------|----------|
| CSD400Y9 | CryoStore 400 Conical EVA Freezing Bag 60-135ml, "Y9" Tube Set | Stock | 12 |
| CSD2500Y9 | CryoStore 2500 Conical EVA Freezing Bag, 730-810ml, "Y9" Tube Set | Non-Stock | 12 |
| CSD2500Y14 | CryoStore 2500 Conical EVA Freezing Bag, 730-810ml, "Y14" Tube Set | Stock | 12 |

*Bag illustrations are not to scale



CSD400Y9

Y9 tube set:
One spike port and one spike with a roller clamp



CSD2500Y14

Y14 tube set:
Two spikes with roller clamps, includes a needle-free valve

OriGen
BIOMEDICAL

7000 Burlleson Rd, Bldg D
Austin, TX USA 78744-3202
Tel: +1 (512) 474 7278
Tel: Toll Free (800) 233 9014
E-mail: sales.us@origen.com