

CRYOSTORE™ Freezing Bags

The OriGen **CRYOSTORE™** bag is intended for freezing blood components.

OriGen
BIOMEDICAL

The industry-preferred freezing bag that has the built-in ease of use, durability, and quality that you need.

EASE OF USE – BUILT-IN

The donor tube is intentionally placed in the corner of the bag for easier debubbling and sealing. OriGen can provide the bag with one of several stock or non-stock tube sets to fit your processing needs. All tubing sets are SCD weldable, allowing for versatile applications. You can order multi-bag sets to reduce processing time and the possibility of contamination in the lab.

NEEDLE-FREE INJECTION PORT

The “N” tubing set has luer-actuated needle-free injection ports. These swabable ports maintain sterility of the sample through numerous actuations (white paper available) and can eliminate the possibility of needle sticks.

ALIQUOTS – MATCHED BAG AND TUBE PRINTING

Each stock bag has a printed number on the label pocket and a donor tube with pre-marked printed segments of approximately 100µL. This can aid in tracing aliquot samples to the bag.

DURABILITY – BETTER MATERIAL, PROVEN PERFORMANCE

The OriGen **CRYOSTORE™** Bag is the industry-preferred cryogenic freezing bag. The hermetically sealed spike ports and film used in every **CRYOSTORE™** bag are made from a robust EVA.

OUR QUALITY STANDARDS

A Certificate of Conformance (CoC) is included with each **CRYOSTORE™** bag. OriGen's **CRYOSTORE™** bags are CE marked and FDA 510(k) cleared. OriGen is certified to ISO 13485:2016, MDSAP certified, and all products are manufactured to GMP guidelines.

ORDERING INFORMATION

To place an online order, please fill out [this form](#).



TECHNICAL INFORMATION

PART CODE	WIDTH (cm*)	LENGTH (cm*)	FREEZE VOLUME MIN-MAX, ML
CS50	7.6	14.7	10-30
CS250	12.7	17.7	30-70
CS500	12.7	22.8	55-100
CS750	12.7	25.7	80-190
CS1000	15.0	28.5	120-275
CS2000	23.0	34.0	180-450

*Length and width are overall bag footprint.

7000 Burleson Rd, Bldg D | Austin, Texas 78744 | 512.474.7278 | Toll Free 800.233.9014 | sales.us@origen.com | origen.com

CE Marked | FDA Cleared | ISO 13485 Certified | MDSAP Certified | Manufactured to GMP Guidelines

PART CODE	TUBE SET DESCRIPTION	PACKAGING	OVERWRAP	DOUBLE POUCHED	CASE QUANTITY	STOCK/ NON-STOCK
CRYOSTORE™ 28 EVA FREEZING BAG, 10-25ML						
CS28V	"V" Tubing Set - Female luer with injection Y to two CS28 bags with pinch clamps	Pouch of 4			24	Non-Stock
CRYOSTORE™ 50 EVA FREEZING BAG, 10-30ML						
CS50	"Standard" Tube Set - Two female luers with caps	Pouch of 4			24	Stock
CS50M3	"M3" Tube Set - One male luer with cap, two female luers with caps, a needle-free injection port and a pinch clamp	Pouch of 4			24	Non-Stock
CS50N	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 4			24	Stock
CS50NS	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 1		✓	12	Stock
CS50S	"Standard" Tube Set - Two female luers with caps	Pouch of 1		✓	12	Stock
CS50V7	Multiple bag set, ten 50 EVA freezing bags, and one 750 EVA freezing bag, "V7" Tube Set - Female luers with caps, "Y" connectors, pinch clamps, slide clamps, co-extruded EVA/PVC tubing and a 30 cm PVC tube	Pouch of 1 Set			2 Sets	Non-Stock
CRYOSTORE™ 250 EVA FREEZING BAG, 30-70ML						
CS250	"Standard" Tube Set - Two female luers with caps	Pouch of 4			24	Stock
CS250F7	"F7" Tube Set - Two female luers with caps and a needle-free valve	Pouch of 4			24	Non-Stock
CS250N	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 4			24	Stock

The OriGen CryoStore™ bag is indicated for freezing blood components. Carefully read Instructions for Use before using this product. Sterile, single use only, made in the USA.

PART CODE	TUBE SET DESCRIPTION	PACKAGING	OVERWRAP	DOUBLE POUCHED	CASE QUANTITY	STOCK/ NON-STOCK
CRYOSTORE™ 250 EVA FREEZING BAG, 30-70ML						
CS250NS	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 1		✓	12	Stock
CS250NSW	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 1	✓	✓	12	Stock
CS250S	"Standard" Tube Set - Two female luers with caps, one bag	Pouch of 1		✓	12	Non-Stock
CRYOSTORE™ 400 EVA FREEZING BAG, 55-100ML, CUSTOM LENGTH 8 INCHES						
CS400	"Standard" Tubing Set - Two female luers with caps	Pouch of 4			24	Non-Stock
CRYOSTORE™ 500 EVA FREEZING BAG, 55-100ML						
CS500	"Standard" Tube Set - Two female luers with caps	Pouch of 4			24	Stock
CS500i2PW	Custom-width 10.9cm bag, "N" Tube set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 1	✓		12	Non-Stock
CS500i3	Custom length 7.5 inches, "Standard" Tubing Set - Two female luers with roller clamps	Pouch of 4			24	Stock
CS500N	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 4			24	Stock
CS500NS	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 1		✓	12	Stock
CS500NSW	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 1	✓	✓	12	Stock
CS500NSWT	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 1	✓	✓	12	Non-Stock

The OriGen CryoStore™ bag is indicated for freezing blood components. Carefully read Instructions for Use before using this product. Sterile, single use only, made in the USA.

PART CODE	TUBE SET DESCRIPTION	PACKAGING	OVERWRAP	DOUBLE POUCHED	CASE QUANTITY	STOCK/ NON-STOCK
CRYOSTORE™ 500 EVA FREEZING BAG, 55-100ML						
CS500S	"Standard" Tube Set - Two female luers with caps	Pouch of 1		✓	12	Non-Stock
CS500V2	Custom 2-Bag set with spikes, "V2" Tube Set - One nylon spike with cover, one male luer with cap, a needle-free injection port and pinch clamps, 2 sets per pouch	Pouch of 2			12	Non-Stock
CS500V8	8-bag set, "V8" Tube Set - One female luer with cap, with "Y" connectors, roller clamps and needle-free "Y"	Pouch of 1 set			3 Sets	Non-Stock
CRYOSTORE™ 750 EVA FREEZING BAG, 80-190ML						
CS750	"Standard" Tube Set - Two female luers with caps	Pouch of 4			24	Stock
CS750N	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 4			24	Stock
CS750NS	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 1		✓	12	Stock
CS750NSW	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 1	✓	✓	12	Stock
CRYOSTORE™ 1000 EVA FREEZING BAG, 120-275ML						
CS1000N	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port with roller clamps	Pouch of 4			24	Stock
CRYOSTORE™ 2000 EVA FREEZING BAG, 180-450ML						
CS2000N	"N" Tube Set - One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 4			24	Stock
CS2000NS	"N" Tube Set One male luer with cap, two female luers with caps, and a needle-free injection port	Pouch of 1		✓	12	Non-Stock

The OriGen CryoStore™ bag is indicated for freezing blood components. Carefully read Instructions for Use before using this product. Sterile, single use only, made in the USA.

ML09.R09: 2022-04

CRYOSTORE™ Freezing Bags

The OriGen CRYOSTORE™ bag is intended for freezing blood components.

OriGen
BIOMEDICAL

EXAMPLE DRAWINGS OF CRYOSTORE™ BAGS AND TUBE SETS




CS28V




CS50



CS50M3



CS50V7



CS250F7



CS250N



CS400



CS250

The OriGen CryoStore™ bag is indicated for freezing blood components. Carefully read Instructions for Use before using this product. Sterile, single use only, made in the USA.

ML60.R00: 2022-07

OriGen
BIOMEDICAL

CRYOSTORE™ Freezing Bags

The OriGen CRYOSTORE™ bag is intended for freezing blood components.

OriGen
BIOMEDICAL

EXAMPLE DRAWINGS OF CRYOSTORE™ BAGS AND TUBE SETS



CS500




CS500N




CS1000N




CS2000N



CS500V2



CS500V8



CS750



CS750N

The OriGen CryoStore™ bag is indicated for freezing blood components. Carefully read Instructions for Use before using this product. Sterile, single use only, made in the USA.

ML60.R00: 2022-07

OriGen
BIOMEDICAL