

CryoPur™ Cryopreservation Solutions



Sterile filtered dimethyl sulfoxide (DMSO) is available in a variety of packages, concentrations, and volumes designed for lab use. 100% and buffered DMSO solutions are packed in polymeric syringes, glass vials or bags. OriGen DMSO products conform to USP and Ph. Eur. requirements.

DMSO and DMSO/Dextran sterile filled syringes are fitted with a stopcock and are shipped with the actuator blocked for safety. They are available in multiple volumes and concentrations for easy, aseptic transfer.

CryoPur™: 100% DMSO

CryoPur™ vials with white caps contain 100% DMSO, are available in 10, 50, 70 and 100ml fill volumes. Caps are flip off/tear off seal.

CryoPur-D™ DMSO/Dextran Solutions

Sterile, filtered DMSO/Dextran solutions are available in a variety of concentrations and formulations. CryoPur-D™ vials with green flip/tear off caps contain 55% w/v DMSO, 5% with Dextran-40, remainder water for injections.

CryoPur™

Cryopreservation Solutions



CryoPur™ : 100% DMSO

Product Code	Description (stock products)	Case Quantity
CP-10	Vial, 10ml filled volume	12
CP-50	Vial, 50ml filled volume	6
CP-70	Vial, 70ml filled volume	6
CP-100	Vial, 100ml filled volume	6
SP-10	10ml in 20ml syringe with stopcock and SCD tubing	1 unit
SP-20	20ml in 20ml syringe with stopcock and SCD tubing	1 unit
SP-30	30ml in 60ml syringe with stopcock and SCD tubing	1 unit
SP-50	50ml in 60ml syringe with stopcock and SCD tubing	1 unit

Ordering Information

CryoPur-D™ : 55% DMSO, 5% Dextran-40 w/v, remainder water for injection

Product Code	Description (stock products)	Case Quantity
CD-07	Vial, 7ml filled volume	12
CD-50	Vial, 50ml filled volume	6
CD-100	Vial, 100ml filled volume	6
SD-07	7.2ml in 20ml syringe with stopcock and SCD tubing	1 unit
SD-10	10ml in 20ml syringe with stopcock and SCD tubing	1 unit

Non-Stock Products Available

Product Code	Description	Case Quantity
SDS-05A-B	5ml 55% DMSO, 5% Dextran 40 w/v, remainder 0.9% saline + 2ml air in 20ml syringe with stopcock and SCD tubing	1 unit
SDS-15A	15ml 55% DMSO, 5% Dextran 40 w/v, remainder 0.9% saline in 20ml syringe with stopcock and SCD tubing	1 unit

DMSO: Store in dark, well ventilated area, freezes (solidifies) at 18 °C. Recommended storage at 20-30°C.

DMSO/Dextran: Store Refrigerated + 2 to + 8°C.

Packaging: 10ml vials 12 to a case, larger vials 6 to a case, syringes sold individually.

Special order volumes and formulations possible. Minimums may apply.

8-12 week lead time on special order items.

For Research use only in the USA.

CryoPur™ is a registered trademark of OriGen Biomedical Inc
© 2014 OriGen Biomedical. All Rights Reserved

OriGen

BIOMEDICAL

Americas & Asia

7000 Burlleson Rd, Bldg D
Austin, TX USA 78744-3202
Tel: +1 (512) 474 7278
Tel: Toll Free (800) 233 9014
Fax: Toll Free (888) 812 8411 or
+1 (512) 617 1503
E-mail: sales.us@origen.com

Europe & Middle East

OriGen Biomedical GmbH
Romerstrasse 14
72393 Burladingen, Germany
Tel: +49-7475-915591
Fax: +49-7475-914718
E-mail: sales.europe@origen.com