

CryoPur™ Cryopreservation Solutions



Sterile filtered dimethyl sulfoxide (DMSO) is available in a variety of packages, concentrations, and volumes designed for lab use. 100% and buffered DMSO solutions are packed in polymeric syringes, glass vials or bags. OriGen DMSO products conform to USP and Ph. Eur. requirements.

DMSO and DMSO/Dextran sterile filled syringes are fitted with a stopcock and are shipped with the actuator blocked for safety. They are available in multiple volumes and concentrations for easy, aseptic transfer.

CryoPur™: 100% DMSO

CryoPur™ vials with white caps contain 100% DMSO, are available in 10, 50, 70 and 100ml fill volumes. Caps are flip off/tear off seal.

CryoPur-D™ DMSO/Dextran Solutions

Sterile, filtered DMSO/Dextran solutions are available in a variety of concentrations and formulations. CryoPur-D™ vials with green flip/tear off caps contain 55% w/v DMSO, 5% with Dextran-40, remainder water for injections.

CryoPur™

Cryopreservation Solutions



Ordering Information

CryoPur™ : 100% DMSO

| Product Code | Description (stock products) | Case Quantity |
|--------------|---|---------------|
| CP-10 | Vial, 10ml filled volume | 12 |
| CP-50 | Vial, 50ml filled volume | 6 |
| CP-70 | Vial, 70ml filled volume | 6 |
| CP-100 | Vial, 100ml filled volume | 6 |
| SP-10 | 10ml in 20ml syringe with stopcock and SCD tubing | 1 unit |
| SP-20 | 20ml in 20ml syringe with stopcock and SCD tubing | 1 unit |
| SP-30 | 30ml in 50ml syringe with stopcock and SCD tubing | 1 unit |
| SP-50 | 50ml in 50ml syringe with stopcock and SCD tubing | 1 unit |

CryoPur-D™ : 55% DMSO, 5% Dextran-40 w/v, remainder water for injection

| Product Code | Description (stock products) | Case Quantity |
|--------------|---|---------------|
| CD-07 | Vial, 7.2ml filled volume | 12 |
| CD-50 | Vial, 50ml filled volume | 6 |
| CD-100 | Vial, 100ml filled volume | 6 |
| SD-10 | 10ml in 20ml syringe with stopcock and SCD tubing | 1 unit |

DMSO: Store in dark, well ventilated area, freezes (solidifies) at 18 °C. Recommended storage at 20-30°C.

DMSO/Dextran: Store Refrigerated + 2 to + 8°C.

Packaging: 10ml vials 12 to a case, larger vials 6 to a case, syringes sold individually.

Special order volumes and formulations possible. Minimums may apply.

For Research use only in the USA.

OriGen
BIOMEDICAL

7000 Burlison Rd, Bldg D
Austin, TX USA 78744-3202
Tel: +1 (512) 474 7278
Tel: Toll Free (800) 233 9014
E-mail: sales.us@origen.com