


Evolve™ Tissue Culture

The OriGen Evolve™ cell culture bags represent the optimal way to harvest, expand and concentrate critical cell cultures. Made of EVO copolymer, the Evolve™ bag contains no plasticizers or chemical residuals which could leach into and alter cell cultures. The larger size bags have hermetically sealed spike ports.

Evolve™ bags are a closed system, which will reduce the possibility of contamination when compared to using a flask or bottle. The Evolve™ bags are gas permeable and provide a high rate of oxygen and carbon dioxide transfer to cells under culture. Additionally water permeability is very low, so water loss is small and humidified incubators are not required.



For availability in your country, please contact your local distributor.

