

SteriZip™ Overpouch

The OriGen SteriZip™ overpouch is intended to provide a safe, clean environment for your CryoStore™ bags when thawing in a water bath.



The applications of a sterile zipper bag are numerous. The SteriZip™ can provide a safe and clean environment for your CryoStore™ bag during thaw and transport throughout your institution. Even use as a closed container to keep your lab supplies.

Each SteriZip™ is inspected, cleaned and pouched in our clean room and gamma sterilized. The bags will arrive sterile and in the zippered position. Simply place your CryoStore™ product inside the SteriZip™ and then re-zip, ensuring you remove enough air so that the SteriZip™ will allow you to massage the CryoStore™ bag while thawing.




SteriZip™

Overpouch

Stock Products

Product Code	Width, (cm)	Height, (cm)	Bags per Pack	For use with the CryoStore™ bag
SZ1823	18	23	10	CS25, CS50
SZ2533	25	33	10	CS250, CS500, CS750
SZ3038	30	38	10	CS1000, CS2000



Additional Options

Product Code	Width, (cm)	Height, (cm)	Bags per Pack/Pouch	Notes
SZ2533-2M	25	33	10	2mil film thickness
SZ2533P-1	25	33	1	Case of 50
SZ2533P-3	25	33	3	Case of 45

Note: These items are non-stock and have a longer lead time when ordering

General Information

The bags are provided sterile and zipped closed until you are ready to use. The bags are 3mil film thickness unless otherwise noted. Custom packaging is available, such as individual aseptic pouches or in quantities other than 10.

Ordering Information:

To order, contact the office nearest you.

OriGen
BIOMEDICAL

7000 Burluson Rd, Bldg D
Austin, TX USA 78744-3202
Tel: +1 (512) 474 7278
Tel: Toll Free (800) 233 9014
E-mail: sales.us@origen.com