

CryoStore™

Conical Freezing Bag

The OriGen CryoStore™ bag is indicated for cryogenically freezing blood components.

Traceability

A permanent lot number is heat welded into the tail of each bag.

Durability

Better Material, Proven Performance

The hermetically sealed spike ports and film used in every CryoStore™ bag are made from a more rugged EVA than that used in other bags. The proof of durability is in our customers' success.

Ease of Use

Designed In

OriGen can provide your choice of several tube sets (see reverse). You can order multi-bags sets prefabricated to reduce processing time and the possibility of contamination in the lab.



Technical Information

Bag Size	Width, (cm)	Length, (cm)	Freeze volume min-max, ml	Fill volume (max), ml
CSD400	12.7	20.3	60-135	420
CSD2500	27.9	33.7	730-810	3000

CryoStore™

Conical Freezing Bag

General Information


Individually packed in Tyvek® pouch, supplied sterile, single use only.

Ordering Information:

For availability in your country, please contact your local distributor.


Product Code	Description	Stock or Non-Stock	Case Qty
CSD400Y9	CryoStore 400 Conical EVA Freezing Bag 60-135ml, "Y9" Tube Set	Stock	12
CSD2500Y9	CryoStore 2500 Conical EVA Freezing Bag, 730-810ml, "Y9" Tube Set	Non-Stock	12
CSD2500Y14	CryoStore 2500 Conical EVA Freezing Bag, 730-810ml, "Y14" Tube Set	Stock	12

*Bag illustrations are not to scale



CSD400Y9

Y9 tube set:
One spike port and one spike with a roller clamp



CSD2500Y14

Y14 tube set:
Two spikes with roller clamps, includes a needle-free valve

OriGen

BIOMEDICAL

Americas & Asia

7000 Burlison Rd, Bldg D
Austin, TX USA 78744-3202
Tel: +1 (512) 474 7278
Tel: Toll Free (800) 233 9014
Fax: Toll Free (888) 812 8411 or
+1 (512) 617 1503
E-mail: sales.us@origen.com

Europe & Middle East

OriGen Biomedical GmbH
Romerstrasse 14
72393 Burladingen, Germany
Tel: +49-7475-915591
Fax: +49-7475-914718
E-mail: sales.europe@origen.com