




Transfer Bags

The OriGen Transfer Bags are intended for the processing of blood and blood components.

The OriGen transfer bags are designed to work for you. Each bag has a range of access points so that you can ensure you can process your products as needed. The multi-bag sets are available for closed cell processing, can be centrifuged and sterile docked. Each bag is labeled for additional security. The bags are made out of PVC in a range of configurations. For availability in your country, please contact your local distributor.



Configuration	Product Code	Description	Case Quantity
	10-200	One 200ml PVC processing bag, with a male luer, a female luer, a needle-free injection port, and a spike port	24
	10-200-2	Two 200ml PVC processing bag, with a male luer, a female luer, a needle-free injection port and a spike port	12
	10-200-CS25	Two 200ml PVC processing bags, with a male luer, a female luer, a needle-free injection port, and a spike port, and a CS25 for final freezing of product	12

Americas & Asia

7000 Burleson Rd, Bldg D
Austin, TX USA 78744-3202
Tel: +1 (512) 474 7278
Tel: Toll Free (800) 233 9014
Fax: Toll Free (888) 812 8411 or
+1 (512) 617 1503
E-mail: sales.us@origen.com

Europe & Middle East

OriGen Biomedical GmbH
Romerstrasse 14
72393 Burladingen, Germany
Tel: +49-7475-915591
Fax: +49-7475-914718
E-mail: sales.europe@origen.com