




# Transfer Bags

The OriGen Transfer Bags are intended for the processing of blood and blood components.

The OriGen transfer bags are designed to work for you. Each bag has a range of access points so that you can process your products as needed. The multi-bag sets are available for closed cell processing, can be centrifuged and sterile docked.\* Each bag is hot stamped for additional traceability. The bags are made of PVC in a range of configurations. For availability in your country, please contact your local distributor.



Configuration	Product Code	Description	Case Quantity
	10-200	One 200ml PVC processing bag, with a male luer, a female luer, a needle-free injection port, and a spike port	24
	10-200-2	Two 200ml PVC processing bag, with a male luer, a female luer, a needle-free injection port and a spike port	12
	10-200-CS25	Two 200ml PVC processing bags, with a male luer, a female luer, a needle-free injection port, and a spike port, and a CS25 for final freezing of product	12

**OriGen**  
BIOMEDICAL

7000 Burleson Rd, Bldg D  
Austin, TX USA 78744-3202  
Tel: +1 (512) 474 7278  
Tel: Toll Free (800) 233 9014  
E-mail: sales.us@origen.com

[www.origen.com](http://www.origen.com)

\* Specific centrifuge parameters have not been validated by OriGen. Customers should validate the bag with their own institution protocols.

© 2014 OriGen Biomedical. All Rights Reserved  
ML24.02.R02; 2019-02