




# PRODUCT CODE DEFINITIONS

# CS250NS



## Product Families\*

<b>CS</b>	CryoStore Freezing Bags
<b>CSD</b>	CryoStore Conical Freezing Bags
<b>CP</b>	CryoPur™ Solutions in Vials
<b>CD</b>	CryoPur™-D Solutions in Vials
<b>SP</b>	CryoPur™ Solutions in Syringes
<b>SD</b>	CryoPur™-D Solutions in Syringes
<b>SDS</b>	CryoPur™-S Solutions in Syringes
<b>PL</b>	PermaLife™ Cell Culture Bags
<b>EV</b>	Evolve™ Cell Culture Bags
<b>OW</b>	O-Wrap™ Bag
<b>TV</b>	TissueVault™ Cell & Tissue Freezing Bags
<b>PFB</b>	Liquivent® Lung Lavage

*\*Excluding accessories*

# Origen<sup>TM</sup>

## BIOMEDICAL

### Freeze Volumes - CryoStore

<b>20, 30</b>	20 - 30ml	<b>07</b>	10ml
<b>25, 28</b>	10 - 25ml	<b>30</b>	40ml
<b>50</b>	10 - 30ml	<b>70</b>	100ml
<b>250</b>	30 - 70ml	<b>120</b>	150ml
<b>500</b>	55 - 100ml	<b>240</b>	300ml
<b>750</b>	80 - 190ml	<b>325</b>	450ml
<b>1000</b>	120 - 275ml	<b>750</b>	1000ml
<b>2000</b>	180 - 450ml	<b>1000</b>	1500ml
		<b>2000</b>	3000ml

### Working Volumes - PermaLife<sup>TM</sup>

### Dimensions - TissueVault<sup>TM</sup>

<b>0909</b>	9x9cm	<b>0909</b>	9x9cm
<b>0918</b>	9x18cm	<b>0918</b>	9x18cm
<b>1430</b>	14x30cm	<b>1430</b>	14x30cm
<b>1950</b>	19x50cm	<b>1936</b>	19x36cm
<b>2436</b>	24x36cm	<b>1950</b>	19x50cm

### Dimensions - O-Wrap<sup>TM</sup>

<b>2436</b>	24x36cm
-------------	---------



## CryoStore Tube Sets

- none** Two female luers with caps
- N** One male luer with cap, two female luers with caps, and a needle-free injection port
- F** Female luer with caps
- Y4** One spike, one male luer with cap and 50cm SCD Tube with a needle-free injection port and pinch clamp
- M3** One male luer with cap, two female luers with caps, a needle-free injection port and a pinch clamp
- V7** Female luers with caps, "Y" connectors, pinch clamps, slide clamps, co-extruded EVA/PVC tubing and a 30 cm PVC tube
- V2** One nylon spike with cover, one male luer with cap, a needle-free injection port and pinch clamps

## CryoStore Packaging

- none** 4 bags per pouch
- P** Single-packed, one pouch
- S** Single-packed, double aseptic pouched
- SW** Single-packed, double aseptic pouched with overwrap

## PermaLife™ Tube Sets

- none** One female luer and one NFV with vent cap
- 2G** One female luer and one NFV
- P** One NFV and one spike port, packed with additional spike port
- P2** One NFV and one 15cm SCD tube
- N or T** Two 15cm SCD tubes