




# PRODUCT CODE DEFINITIONS

# CS250NS



## Product Families\*

<b>CS</b>	CryoStore Freezing Bags (Including FLEX and Multi-Chamber)
<b>CSD</b>	CryoStore Conical Freezing Bags
<b>CP</b>	CryoPur™ Solutions in Vials
<b>CD</b>	CryoPur™-D Solutions in Vials
<b>SP</b>	CryoPur™ Solutions in Syringes
<b>SD</b>	CryoPur™-D Solutions in Syringes
<b>SDS</b>	CryoPur™-S Solutions in Syringes
<b>PL</b>	PermaLife™ Cell Culture Bags
<b>EV</b>	Evolve Cell Culture Bags
<b>OW</b>	O-Wrap™ Bag
<b>TV</b>	TissueVault™ Cell & Tissue Freezing Bags
<b>PFB</b>	Liquivent® Lung Lavage

*\*Excluding accessories*



**Freeze Volumes**  
- **CryoStore Bags**

<b>20, 30</b>	20 - 30ml	<b>07</b>	10ml
<b>25, 28</b>	10 - 25ml	<b>30</b>	40ml
<b>50</b>	10 - 30ml	<b>70</b>	100ml
<b>250</b>	30 - 70ml	<b>120</b>	150ml
<b>500</b>	55 - 100ml	<b>240</b>	300ml
<b>750</b>	80 - 190ml	<b>325</b>	450ml
<b>1000</b>	120 - 275ml	<b>750</b>	1000ml
<b>2000</b>	180 - 450ml	<b>1000</b>	1500ml
		<b>2000</b>	3000ml

**Working Volumes**  
- **PermaLife™ Bags**

**Dimensions**  
- **TissueVault™ Bags**

<b>0909</b>	9x9cm
<b>0918</b>	9x18cm
<b>1430</b>	14x30cm
<b>1950</b>	19x50cm
<b>2436</b>	24x36cm

**Dimensions**  
- **O-Wrap™ Bags**

<b>0909</b>	9x9cm
<b>0918</b>	9x18cm
<b>1430</b>	14x30cm
<b>1936</b>	19x36cm
<b>1950</b>	19x50cm
<b>2436</b>	24x36cm



## CryoStore Tube Sets

<b>none</b>	Two female luers with caps; two roller clamps
<b>N</b>	One male luer with cap; two female luers with caps and three roller clamps; a needle-free injection port
<b>F</b>	Female luer with cap
<b>Y4</b>	One spike with cover; one male luer with cap and 50cm SCD tubing with a needle-free injection port and pinch clamp
<b>M3</b>	One male luer with cap; two female luers with caps and three roller clamps; a needle-free injection port and a pinch clamp
<b>V7</b>	Multi-bag Set: Female luer with cap; "Y" connectors; pinch clamp; slide clamps; co-extruded EVA/PVC tubing and a 30 cm PVC tube
<b>V2</b>	Double Bag Assembly: One spike with cover; one male luer with cap; a needle-free injection port and three pinch clamps

## CryoStore Packaging

<b>none</b>	4 bags per pouch
<b>P</b>	Single-packed, one pouch
<b>S</b>	Single-packed, double aseptic pouched
<b>SW</b>	Single-packed, double aseptic pouched with overwrap

## PermaLife™ Tube Sets

<b>none</b>	One female luer and one NFV with vent cap
<b>-2G</b>	One female luer and one NFV
<b>P</b>	One NFV and one spike port, packed with additional spike port
<b>P2</b>	One NFV and one 30cm SCD tube packed with female luer cap
<b>T</b>	Two 15cm SCD tubes, one sealed and one with porous luer plug