




PRODUCT CODE DEFINITIONS

CS250NS



Product Families*

CS	CryoStore Freezing Bags (Including FLEX and Multi-Chamber)
CSD	CryoStore Conical Freezing Bags
CP	CryoPur® Solutions in Vials
CD	CryoPur-D Solutions in Vials
SP	CryoPur Solutions in Syringes
SD	CryoPur-D Solutions in Syringes
SDS	CryoPur-S Solutions in Syringes
PL	PermaLife® Cell Culture Bags
EV	Evolve Cell Culture Bags
OW	O-Wrap® Bag
TV	TissueVault™ Cell & Tissue Freezing Bags
PFB	Liquivent® Lung Lavage

**Excluding accessories*



Freeze Volumes - CryoStore Bags

20, 30	20 - 30ml	07	10ml
25, 28	10 - 25ml	30	40ml
50	10 - 30ml	70	100ml
250	30 - 70ml	120	150ml
500	55 - 100ml	240	300ml
750	80 - 190ml	325	450ml
1000	120 - 275ml	750	1000ml
2000	180 - 450ml	1000	1500ml
		2000	3000ml

Working Volumes - PermaLife® Bags

Dimensions - TissueVault™ Bags

0909	9x9cm	0909	9x9cm
0918	9x18cm	0918	9x18cm
1430	14x30cm	1430	14x30cm
1950	19x50cm	1936	19x36cm
2436	24x36cm	1950	19x50cm

Dimensions - O-Wrap® Bags

0909	9x9cm
0918	9x18cm
1430	14x30cm
1936	19x36cm
1950	19x50cm
2436	24x36cm



CryoStore Tube Sets

none	Two female luers with caps; two roller clamps
N	One male luer with cap; two female luers with caps and three roller clamps; a needle-free injection port
F	Female luer with cap
Y4	One spike with cover; one male luer with cap and 50cm SCD tubing with a needle-free injection port and pinch clamp
M3	One male luer with cap; two female luers with caps and three roller clamps; a needle-free injection port and a pinch clamp
V7	Multi-bag Set: Female luer with cap; "Y" connectors; pinch clamp; slide clamps; co-extruded EVA/PVC tubing and a 30 cm PVC tube
V2	Double Bag Assembly: One spike with cover; one male luer with cap; a needle-free injection port and three pinch clamps

CryoStore Packaging

none	4 bags per pouch
P	Single-packed, one pouch
S	Single-packed, double aseptic pouched
SW	Single-packed, double aseptic pouched with overwrap

PermaLife® Tube Sets

none	One female luer and one NFV with vent cap
-2G	One female luer and one NFV
P	One NFV and one spike port, packed with additional spike port
P2	One NFV and one 30cm SCD tube packed with female luer cap
T	Two 15cm SCD tubes, one sealed and one with porous luer plug